



Description for Microcentrifuge:

From the intuitive controls, to the innovative Thermo Scientific ClickSeal™ biocontainment lid, the Sorvall Legend Micro 21 microcentrifuge delivers the perfect blend of capabilities to support your micro-volume protocols such as nucleic acid or protein lysate preparation and PCR reaction set-up—all in a remarkably small footprint.

The Sorvall Legend 21 microcentrifuge provides a new separation experience with:

Exceptional acceleration up to 21,000 x g in 12 seconds

Refrigerated versions offer cooling in just 9 minutes

Holds 36 x 0.5mL microtubes, 24 x 2mL tubes or 8 x 8 PCR strips

ClickSeal biocontainment rotor lid “clicks” in place, ensuring securely contained samples

New PPSU lids have improved strength and chemical resistance

Intuitive controls and easy to read display simplifies operation

Sleek, space-saving footprint

Rugged design for years of dependable service

Quiet performance provides safe, stress free environment

Description for Corning™ LSE™ Vortex Mixer:

Unit features a powerful, counter-balanced motor drive, continuous or “touch” mix modes. A range of accessories are available for all common sample tubes and containers. Variable speed for gentle to vigorous mixing. Optimized counter balance for minimal vibration. Lightweight and portable. Mixes a variety of tubes and plates with optional accessory heads. Touch and continuous mix modes.